



REANIMAT[®]-O₂-Blend

Insufflation and Suction

REANIMAT[®]-O₂-Blend


Easy, safe and reliable

The supply of oxygen-air blends by so called blenders are part of the clinical standard today. They are often combined with a connected flowmeter for dosing the flow.

The insufflation assembly REANIMAT-O₂-Blend is a universal device for the supply of medical oxygen-air-blends both for stationary and mobile use.

Dosing of the flow is possible from 0 to 60 l/min, according to the intended clinical use.

Optionally a suction unit can be integrated in the insufflation assembly which makes this device indispensable for various non-invasive applications of breathing gas support in pediatrics, neonatology and adult care.



Wide range of application

Due to its modular configuration the assembly can be used in all fields where oxygenated medical blended gas is required. A combination with non-invasive respiration equipment is usually possible, provided the respiration equipment has a pressure limiting facility.

Application examples: Oxygenation of respiration units such as

- Resuscitation bag
- Emergency respirators
- CPAP devices

Second take off as option

It is an ideal possibility to have available oxygenated blended gas from a second take off.

By using a Y-piece, the sum of both flows can be made available. In this way it is possible to use one take-off for the rough selection and the second take off for the fine tuning.

Example: 0-15 l/min for rough selection and 0-6 l/min for fine tuning, or, 0-30 l/min for rough selection and 0-15 l/min for fine tuning.



Avoiding needless gas consumption

Adequate devices being offered in the market usually consume 8-13 l/min of blended gas, even if the flow valves are closed and the device is obviously out of service. Such consumption cause costs and is avoidable.

The insufflation assembly REANIMAT-O₂-Blend has an on/off switch by which the gas supply can be interrupted after use. It is not necessary to put the wall plugs into park position or to close the cylinder valves. It is sufficient to simply switch off the REANIMAT-O₂-Blend.

Hence the device will pay for itself within short time.



Selected parameters do not change

The selected parameters are of impressive consistency. If the flow is changed, the oxygen concentration remains unchanged, and if the oxygen concentration is changed, the selected flow remains constant.

The flow and concentration values remain unchanged after switching off the device. It is not necessary to close the flow valves after use – the selected values are available again when switching on the device next time.




Integrated suction as option

In general it is necessary to aspirate the airway in regular intervals. Due to its modular configuration this insufflation assembly is optionally available with an integrated suction unit.


Same as the insufflation assembly, also the suction unit has an on/off switch. So the selected vacuum remains constant when switching off the device. When switching on the device next time, it is not necessary to adjust the vacuum again.

By this convenient device combination of insufflation and suction a compact device is available which has only one each pressure gas connection for oxygen and Air. Gas sockets are not blocked unnecessarily which also contributes to a clear arrangement of the patient's care place.


Various connection alternatives



Air and oxygen from central gas pipeline



Supply from central gas and cylinder




Supply from two cylinders air and oxygen

Automatic changeover to emergency supply possible:

In many cases it is essential to secure the gas supply for every situation.

By interconnecting the preference valves WY1475 and WY1476, that are available as accessories, the pressure gas supply through gas cylinders is automatically secured, in case the permanent compressed gas supply fails or is interrupted, for example when the patient has to be moved from one place to another place within the hospital.



REANIMAT[®]-O₂-Blend

Ordering information

Insufflation assembly 0-6 l/min with blender 21-100 % oxygen	WY1480
Insufflation assembly 0-15 l/min with blender 21-100 % oxygen	WY1481
Insufflation assembly 0-30 l/min with blender 21-100 % oxygen	WY1482
2nd take-off 0-6 l/min for WY1480 / WY1481 / WY1482	WY1485
2nd take-off 0-15 l/min for WY1480 / WY1481 / WY1482	WY1486
2nd take-off 0-30 l/min for WY1480 / WY1481 / WY1482	WY1487
Suction unit 0-0.9 bar vacuum	WY1490

Accessories

Oxygen preference valve, 2 inputs - 1 take-off	WY1475
Air preference valve, 2 inputs - 1 take-off	WY1476
Y-piece for inter-connecting two blended gas take-offs	WY1478

Connection hoses as per ISO 13260 with gas specific NIST unions are included in delivery as standard. Country-specific connection hoses are available on request.

Technical data

Pressure gas supply	Medical oxygen, 3.4 to 5 bar (50-73 psi) Pure medical compressed air, 3.4 to 5 bar (50-73 psi)
Connection thread for pressure gas	Gas-specifically NIST according to DIN EN 739
Blender	21% to 100% ± 3% Fi=2 accuracy
Gas failure alarm	1.5 to 2 bar between the inlet pressures
Gas consumption oxygen and air	1 to 8 l/min. plus selected take-off
Gas consumption for the suction unit	0 to 30 l/min
Vacuum of suction unit	0 to -0.9 bar vacuum, infinitely variable
Blended gas take-off base device	WY1480: 0-6 l/min WY1481: 0-15 l/min WY1482: 0-30 l/min
Second take of blended gas (option)	WY1485: 0-6 l/min WY1486: 0-15 l/min WY1487: 0-30 l/min
Classification according Directive 93/42/EWG	Ila
Certification	CE 0197



Weyer GmbH

Herrenhöhe 4, D-51515 Kürten-Herweg / Germany

Tel.: (+49) 02207 9698-0

Fax.: (+49) 02207 4750

E-Mail: info@weyermed.com

www.weyermed.com