

Baby changing and examination



## **THERMOCARE** *Convenience*

Table top heating module providing warmth and comfort

# THERMOCARE *Convenience*

New born babies require a warm and cosy environment. Especially preterm babies who have an extended and increased need of devotion, warmth and shelter. Therefore the examination and treatment area should be furnished for these small patients to feel as comfortable as possible. The baby warmer of our successful product line THERMOCARE helps to achieve this goal. The soft and warm gel pad reminds the baby of the contact with the parents' skin and has a calming and comforting effect.

The heating module's compact design allows placement on any suitable table top and an easy integration into baby changing areas without major installation.

Combined with a CERAMOTHERM® radiant warmer, integral workstations will provide ambient conditions equivalent to a mobile resuscitation unit.

The comfortable conductive under-pad heating warms up even hypothermic infants quickly and efficiently. The temperature is selected in increments of 0.1 °C, whilst displayed and monitored electronically. Alarms will alert in case of high or low temperature, system faults or power failure. A separate safety circuit will cut off the heating in the case of a fault.

## Ordering information

THERMOCARE Convenience,  
Baby warmer with gel pad,  
4 rubber feet and integrated storage place  
Order No. WY2011



## Technical Data

Dimensions	715 x 600 x 40 resp. 100 mm (D x W x H)
Weight	12 kg
Operating voltage / power supply	230 V~ / 50 60 Hz
Max. power input	160 W / 0.7 A
Temperature selection range, standard	35 °C to 37 °C
Temperature selection range, extended	30 °C to 38.5 °C
Increments of temperature selection	0.1 °C
Resolution of temperature selection	0.1 °C
Display of temperature selection	30 °C to 38.5 °C
Display of actual temperature	0 °C to 50 °C
Display of residual heat	Down to 30 °C, then automatic switch-off
Safety cut-off	At 40.5 °C
Protection class	1
Type complete device	B
Type under-pad heating	BF
MDD-class	Ila



### Weyer GmbH

Herrenhöhe 4, D-51515 Kürten-Herweg  
Tel.: (+49) 02207 9698-0  
Fax.: (+49) 02207 4750  
E-Mail: [info@weyermed.com](mailto:info@weyermed.com)  
[www.weyermed.com](http://www.weyermed.com)

